

# Effect of gaze direction and lid squeezing on intraocular pressure in glaucoma patients using a novel implantable intraocular sensor

J. van den Bosch<sup>1,2</sup>, V. Pennisi<sup>2</sup>, H. Thieme<sup>2</sup>, L. Choritz<sup>2</sup>

<sup>1</sup>Implandata Ophthalmic Products GmbH, Hannover, Germany

<sup>2</sup>University Eye Clinic, Otto-von-Guericke University Magdeburg, Germany

CONTACT J. van den Bosch

[jvandenbosch@implandata.com](mailto:jvandenbosch@implandata.com)

## Purpose


Quantification of immediate intraocular pressure ('IOP') changes in response to voluntary eye globe movement and lid squeezing

## Methods

- Eight glaucoma patients implanted with the EYEMATE-IO IOP sensor, measured every 0.1 second via an external antenna around the study eye (Fig. 1).
- Seated towards a static grid with determined gaze angle positions (Fig. 2).
- Clockwise directed gaze positions in 8 directions, intervals alternated by primary position as baseline measurement. Circles completed in increased eccentricity (10, 20 and 25° , Fig. 3).
- IOP fluctuations depicted as mean IOP change (Fig. 4-6).
- Instructed lid squeezing alternated by primary position, depicted as peak IOP change (Fig. 7, 8).




1 EYEMATE™ sensor. Implantation during cataract surgery in the ciliary sulcus, in front of the intraocular lens.



2 Set-up with the patient in a chinrest 2.5 meters in front of the Harms wall, fixating at the center as primary position.


## Results

### Continuous IOP recordings during different gaze positions




3 Example of 8 different gaze positions in 12-second intervals following a 25° circle. Primary position is indicated with a dot. Upward gazes provoke IOP increases up to 10 mmHg.

### Analysis




4 Change in IOP expressed as the mean IOP during the gaze position compared to the mean of the previous primary position.

### Intersubject variability in mean IOP changes




5 Two individual heat maps (study eye is the right eye). Subject specific differences in IOP range and changes in respect to gaze positions are present (color legend Fig. 6).

### Combined overview per gaze position




6 Heat map of the combined change in mean IOP per gaze direction (mean of 3 repetitions). Vertical and temporal gazes increase IOP while nasal and downward gazes decrease IOP.

### IOP recordings of lid squeezing



7 Example lid squeeze, alternated by 10-second baseline intervals. The ocular pulse amplitude is apparent. Change in IOP expressed as the peak IOP during lid squeezing compared to the previous baseline mean.

### Intersubject variability in peak IOP change



8 Peak IOP change during lid squeezing in 8 patients, averaged over 5 repetitions (mean as green circle; lines as minimum and maximum values).

## Conclusions

- The Eyemate-IO IOP sensor enables short term observation of IOP fluctuations in glaucoma patients without movement restriction.
- Asymmetry in IOP response in different gaze directions and profound IOP increases during lid squeezing, in line with previous literature.<sup>1,2</sup> Extend of observed intersubject variability has not been described.
- Future analyses should explore ocular pulse amplitude as well as aqueous humor outflow dynamics in response to controlled increases in IOP.

## References

1. Nadine, M. Intraocular pressure changes in secondary positions of gaze in normal subjects and in restrictive ocular motility disorders. Graefes Arch Clin Exp Ophthalmol. 1988; 226: 8-10.
2. Coleman, J. Trokel, S. Direct-recorded Intraocular Pressure Variations in a Human Subject. Arch Ophthalmol. 1969; 82: 637-641.